

Triple-MRE-seq library construction protocol v5 Costello Lab March 2011

Notes:

- A. For all Qiagen gel extraction steps (Qiaquick and MinElute), melt gel slice at 37° C instead of 50° (see Quail et al 2008 Nature Methods). Incubate for 30 min, vortexing frequently - make sure gel is completely melted before proceeding.
- B. For Qiagen Qiaquick and MinElute column purifications (not gel extractions) use Qiagen buffer PB, not PBI (newer Qiagen kits don't even have PBI).
- C. For size-selection after adapter ligation, be aware that the adapters are only partially double-stranded so adapter-ligated fragments might run at slightly different position than expected.
- D. All gel purifications are performed with TAE agarose gels. We add 1:10000 Sybr-Safe to gel after microwaving (instead of Et Br).
- E. The Illumina Genomic DNA Sample kit manual specifies Bio-Rad certified low-range Ultra Agarose for gel extraction steps; the ChIP-seq Sample kit manual does not specify a particular agarose. I have been using the Bio-Rad Ultra Agarose.
- F. Minimum amount of starting DNA is currently 1.5 µg.

PROTOCOL:

1. Start with high-quality, high molecular weight gDNA. It is very important to prevent random shearing of DNA. Isolate DNA using the Costello lab's "genomic DNA extraction from frozen tissue" protocol. Include overnight proteinase K digestion and a **1 hour**, 37 degree DNase-free RNase treatment (Roche catalog# 11 119 915 001) (final RNase concentration = 40 µg/ml). Perform 2 PCI and 2 chloroform extractions using 2 ml "Light" phase lock gels. Be careful at all steps not to shear DNA – always use wide-bore tips and pipet gently. Resuspend DNA in TE (not water). Store DNA at 4° C.
2. Digest DNA using Fermentas methylation-sensitive enzymes *HpaII*, *Hin6I*, *AciI*. Each digest is performed separately; use 0.5 µg DNA per digest. Use 10 units of enzyme per µg of DNA. Add half of enzyme, incubate 3 hrs at 37 ° C, then add rest of enzyme, incubate additional 3 hrs at 37° C. Total digest time = 6 hrs.

Setup for each digest (0.5 µg DNA each : USE WIDE-BORE FILTERED TIPS)

| | |
|---------------------------------------|--|
| 10X buffer | 5 µl |
| Fermentas restriction enzyme, 10 U/µl | 0.25 µl, then additional 0.25 µl after 3 hrs |
| DNA in TE | add volume for 0.5 µg |
| Illumina ultrapure water | to 50 µl |

mix reaction gently – stir, tap gently, or pipet gently with wide-bore tips

3. Clean up digest and resuspend in smaller volume. (Use **wide-bore filtered tips**). Add 1 volume of phenol/chloroform/isoamyl alcohol to each digest and mix. Next, combine digests into 1 phase lock gel tube (MREs should be inactivated by addition of PCI). Continue with PCI extraction followed by one chloroform extraction, using phase lock gel. Precipitate DNA by salt/ethanol precipitation (1/10 volume of 3M sodium acetate, pH 5.2; 2.5X volume of 100% EtOH) with 1 µl glycogen. After addition of EtOH, place in -20 overnight. Resuspend pellet in 12 µl Qiagen EB buffer. Nanodrop 2 µl to check concentration and purity. Save remainder of digests at -20.

4. Make a 50 ml, 2.5% TAE gel (Owl small gel; use thick 10-well comb). **Using wide-bore filtered tips** mix DNA with 4X tris/EDTA/sucrose loading buffer (does not contain dyes). Loading gel: skip first and last wells. Load combined digests on one side of gel; load 5 μ l of 100 bp ladder and 3 μ l of NEB Low MW ladder on other side. Load a lane next to marker with 0.5 μ l dye only (BPB and XC) to help monitor gel running. Run at 120 V 1 hr or longer, depending on size-selection.

5. Size selection of digest: Cut out 50-300 bp (and additional size fractions if wanted) with clean scalpel. Be careful to exclude fragments <50 bp. Minimize exposure to UV, and take a picture before and after excising gel slice. Try to make cut as horizontal as possible, and minimize the size of each gel slice. Purify each gel slice with a single Qiagen Qiaquick column (for gel slices < 400 mg) or two Qiagen MinElute columns (for gel slices > 400 mg). Melt gel at 37° C for ~ 15 min with frequent mixing. Elute DNA in a total of 32 μ l of Qiagen EB buffer. Concentration will be low. After this step, wide-bore tips are not necessary.

6. Repair reaction. Use all 30 μ l of size-selected digest. In contrast to standard Illumina Genomic DNA kit protocol, **exclude** T4 DNA polymerase and T4 PNK. Enzymes and buffers are from Illumina Genomic DNA kit. Make Klenow DNA polymerase dilution: 1 μ l + 4 μ l Illumina water. Discard unused diluted Klenow after use.

| | |
|--|------------|
| Illumina water | 12 μ l |
| DNA sample | 30 |
| T4 DNA ligase buffer with 10mM ATP | 5 |
| dNTPs | 2 |
| Klenow DNA polymerase, diluted 1:5 in H ₂ O | <u>1</u> |
| total volume | 50 |

incubate in thermal cycler 30 min at 20°C

7. Clean up reaction with Qiagen **Qiaquick** column, as per Qiagen instructions. Elute in 32 μ l EB.

8. Addition of 3' A (using Illumina reagents)

| | |
|------------------------|------------|
| DNA sample | 32 μ l |
| Klenow buffer | 5 |
| dATP | 10 |
| Klenow 3'-5' exo minus | <u>3</u> |
| total | 50 |

incubate in thermal cycler 30 min at 37°C

9. Clean up with Qiagen **MinElute** column, elute in 10 μ l EB.

10. Adapter Ligation. Use 1 μ l of 1:10 dilution of PE adapters. Dilute 1 μ l Illumina adapter oligo mix in 9 μ l Illumina Ultrapure water. Store extra diluted adapter oligo mix at -20 °C.

| | |
|-------------------|------------|
| DNA sample | 10 μ l |
| DNA ligase buffer | 15 |

| | |
|--|-----------|
| 1:10 diluted PE adapter oligo mix | 1 |
| DNA ligase | <u>4</u> |
| total | 30 |

incubate 15 min at room temperature

11. Clean up with Qiagen **MinElute** column, elute in 20 µl EB.

12. Size-select adapter-ligated fragments. Follow Illumina Chip-seq kit protocol for making and loading 50 ml 1X TAE gel (Owl small gel; use thick 10-well comb).

Mix 20 µl DNA with 6 µl 4X tris/EDTA/sucrose loading buffer (does not contain dyes). Loading gel: skip first and last wells. Load all 26 µl sample on one side of gel; load 5 µl of 100 bp ladder and 3 µl of NEB Low MW ladder on other side of gel. Load a lane next to markers with 0.5 µl dye only (BPB and XC) to help monitor gel running. Run at 120 V 1 hr. Cut out gel slice corresponding to original size selection + adaptors. Optional: excise an additional gel slice of the same size from a lane without sample; use this as negative control for PCR. Gel-purify each slice with a single Qiagen **Qiaquick** column (for gel slices < 400 mg) or two or more Qiagen **MinElute** columns (for gel slices > 400 mg). Elute in 30 µl EB. Use 1/3 (10 µl) for PCR; save rest at -20 °C.

13. PCR enrichment of adapter-modified fragments. Can use water and/or blank gel slice as negative control.

| | |
|---------------------------------|-----------|
| DNA | 10 |
| Illumina Phusion DNA polymerase | 25 |
| Illumina PCR primer 1.1 | 1 |
| Illumina PCR primer 2.1 | 1 |
| Illumina Ultrapure water | <u>13</u> |
| total | 50 |

cycling conditions:

30 sec 98 °C

15 cycles of:

10 sec 98 °C

30 sec 65 °C

30 sec 72 °C

Then: 5 min 72 °C

Then: hold at 4 °C

14. Size select 170-420 bp by Qiagen **MinElute** gel purification. Elute in 16 µl EB.

→ 2 µl for nanodrop

→ 1 µl for Bioanalyzer

→ ? µl for subcloning and Sanger sequencing

15. If subcloning: use Invitrogen pCRII-Blunt Topo cloning kit. Try for 3:1 ratio of insert:vector for ligation. Sequence ~ 25 colonies per sample.